Medtrum EasyView Website

Korisničko uputstvo

Za negovatelje

Humalog® je trgovačka marka Eli Lilly i Kompanije.

NovoRapid® je trgovačka marka Novo Nordisk A/S.

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


Nakon što se prijavite, videće Vašu listu pacijenata.

Name	SG (mmol/L)	Basal (iU/L)	Active Insulin (iU)	Action
Alan				Stop following
Yves				Stop following
Luc	8.3		0	Stop following
Mal	6.9	0.8	0	Stop following
Dele	8.1	0.9	0	Stop following
Clair				Stop following

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
Connections		Add Connection		
Name	SG (g/lmmol/L)	Basal (U/h)	Active Insulin (U)	Action
Alan				Stop following
Voss				Stop following
Eva	8.3		0	Stop following
Nell	6.0	0.8	0	Stop following
Carla	8.1	0.9	0	Stop following
Clark				Stop following

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A form titled 'Patient profile' with a light blue header. The form contains the following fields: 'Username: test18@medtrum.cn', 'Photo: [person icon]', 'Full Name: Clair Simpson', 'Alias: [text input]', 'Mobile: --', 'Telephone: --', and 'Email: --'. At the bottom right, there are two buttons: 'Save' and 'Cancel'.

Kliknite na pacijenta koga želite da pogledate u listi pacijenata, i tako ćete prikazati njen/njegov ekran za praćenje.

Kliknite na **Settings (Podešavanja)** na plavoj liniji kako biste prikazali ekran podešavanja.

Report Settings
Save Cancel

Report Settings (mmol/L)

Show alerts On Off

Glucose Target Range

Overall 4.5 - 8.8

Hypoglycemia 3.1

		Glucose Target Range		Time Period		Post Meal Analysis	
		Low	High	From	To	From	To
Breakfast	Before	4.6	6.5	06:15	07:00	01:00	01:30
	After	5.8	8.7				
Lunch	Before	4.6	6.5	11:00	12:00	01:00	01:30
	After	5.8	8.7				
Dinner	Before	4.6	6.5	16:00	17:00	01:00	01:30
	After	5.8	8.7				
Night-time		4.6	6.5	21:00	05:30		

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



When you change the end date on one of the reports, the end dates on all other reports will be synchronized automatically.

7.1 Dashboard

The dashboard report is used for reviewing the statistics of CGM readings and insulin pump delivery. It provides detailed figures for analyzing and managing diabetes.


The first part of the report shows the overall glucose statistics: average BG, CGM reading, daily insulin delivery, carbs-intake, exercise time as well as standard deviation (SD) and number of input (#). An estimated A1C is calculated from SG values.

Glucose Statistics										
BG		CGM		Insulin		Carbs		Exercise		
Average		Average		Average Daily Dose		Average Carbs/Day		Average time/Day		
9.4 mmol/L		9.1 mmol/L		37.33 U		44 g		---		
SD	#	SD	#	SD	# days	SD	#	SD	#	#
3.17	30	3.66	4901	2.57	4	18.6	6	--	--	--


Estimated A1C: 7.3% calculated from SG values Note: Estimated A1C does not replace lab measurement

The second part displays the percentile report, in which you can view the SG distribution at certain daily time points during the selected days.


Reports

The third part of the report includes detailed glucose and Insulin delivery statistics. You can also view the glucose and insulin distribution of the selected days in the pie chart.

You can set the period by clicking on **Set**, and then selecting the duration and the end date.

A dialog box titled "Set period" with a close button (X) in the top right corner. It contains a "Duration" dropdown menu set to "7 days", a "From" date field set to "19-10-2017", and a "To" date field set to "25-10-2017" with a calendar icon. At the bottom, there are two buttons: "Apply" and "Cancel".

You can print this report by simply clicking on **Print**.

7.2 Daily summary

This report presents the sensor glucose curve, calibration marks, basal rates, bolus doses, alerts, carbohydrate records, and exercise records of a selected day.

Reports

By moving the mouse pointer along the curves, you can see more detailed information.

The default date is always the last date before the day you open the website.

Date < 24/10 2017 >

You can change the date by clicking on Date, and choose the date on the calendar.

You can print this report by simply clicking on **Print**.

7.3 Day by day overview


This report presents the meter blood glucose readings, sensor glucose curves, insulin delivery curves, and carbs intake of 7 or 14 days. The default period is 7 days.

By moving the mouse pointer along the curves, you can see more detailed information.

Reports

You can set the period by clicking on **Set**, and then selecting the duration (7 or 14 days) and the end date.

Set period


Duration 7 days

From 19-10-2017 To 25-10-2017


Apply Cancel

You can print this report by simply clicking on **Print**.


7.4 Sensor overlay

This report displays the sensor data overlay and meter blood glucose readings within a given number of days (maximum 7 days) before a selected date. All the daily SG curves and meter blood glucose readings are displayed in an overlap graph so that you can easily see the pattern of glucose levels in a given period. If you print it, you will also see the daily average SG, maximum SG, minimum SG, and the number of hypoglycemia episodes.

By moving the mouse pointer along the curves, you can see more detailed information.

You can set the period by clicking on **Set**, and then selecting the duration, the date scale, the end date and the days you want to see.

You can print this report by simply clicking on **Print**.

7.5 Trend analysis

This report shows the trend graph and statistics within a given number of days (maximum 90 days).

The Glucose Trend graph and statistics include:

- Basic statistics
 - Avg.: The average of all the sensor readings.
 - Max: The highest sensor reading.
 - Min: The lowest sensor reading.
 - Med: The Midpoint of all the sensor readings.
 - Q1: 25% of sensor readings are lower than this value.
 - Q3: 75% of sensor readings are lower than this value.
 - IQR: The interquartile range (IQR) is the difference between the Q1 and Q3.
 - SD: The Standard deviation (SD) measures the amount of variability or dispersion, from the individual sensor readings to the average.
 - CV: The Coefficient of variation (CV) is a measure of relatively variability.
 - SE Mean: The Standard Error (SE) Mean measures how spread out the sensor readings are.
- Distribution %: The distribution of the sensor readings. The high and low limits are set in the report settings.
- Stability:
 - AUC above limit: The Area under a Curve (AUC) above limit is a measure of the high glucose above the high limit.

- AUC below limit: The AUC below limit is a measure of the low glucose below the low limit.
 - # of high / low: Shows how many times the sensor readings are higher than the high limit or lower than the low limit.
- BG statistics
- Avg: The average of all the BG readings.
 - SD: The Standard deviation (SD) measures the amount of variability or dispersion, from the individual BG reading to the average.
 - # of Hypo: Shows how many times the BG readings are lower than 3.1 mmol/L (56 mg/dL).
 - BG readings: The number of the BG readings.

The Insulin trend graph and statistics include:

- Avg. Daily Total: The average of total insulin delivered per day.
- Avg. Daily Basal: The average of daily Basal delivered per day.
- Avg. Daily Bolus: The average of daily Bolus delivered per day.
- Avg. # of bolus/Day: The average number of times to deliver a bolus per day.
- Avg. bolus(U)/Each: The average insulin delivered for each Bolus.

The Carbs statistics include:

- Carbs (g): The carbohydrates intake
- # of meals: The number of meals.

Reports

The Exercise statistics include:

- Duration (min): The total duration of exercise
- # of Exercises: The number of exercises .


By moving the mouse pointer on the graphs, you can see more detailed information.

You can change the period by clicking on **Set**, and then selecting the start date, end date and the days you want to see.

You can print this report by simply clicking on **Print**.


7.6 Events

You can view here both the events recorded on patient's EasyTouch/EasySense app and those recorded in patient's PDM if the patient has uploaded PDM data to the cloud.


				Period: 12/11/2019 - 18/11/2019, All		Set	Print
	Time	Event Type	Details	Notes	Source	Action	
<input type="checkbox"/>	18-11-2019 21:29:00	BG	12.8 mmol/L		PDM	Edit	Delete
<input type="checkbox"/>	18-11-2019 21:29:00	Carbs	66 g		PDM	Edit	Delete
<input type="checkbox"/>	18-11-2019 21:29:00	Exercise	25 min Light		PDM	Edit	Delete
<input type="checkbox"/>	18-11-2019 13:49:00	Injection	10 U		app	Edit	Delete
<input type="checkbox"/>	18-11-2019 13:49:00	Injection	10 U		app	Edit	Delete

To create a new event from the website, click on **Create**, record the event in the pop-up, and then click on **Save**.


To edit an existing event, click on **Edit**, make changes in the pop-up, and then click on **Save**.

Reports


To delete events from the website, tick them and then click on **Delete**.


You can set the period by clicking on **Set**, and then selecting the end date and the number of days you want to see.


If you want to view only one category of events, click on **Set** and use the event filter.

Click on **Quick Print** in the blue bar to go to the Quick Print screen.

On this screen you can quickly print different reports: select the duration, the start date, the end date, and the days you want to see, tick the reports you want to print, and click on **Create Report** to print all of them.


The screenshot shows the 'Quick Print' configuration interface. It includes a 'Quick Print Settings' section with a 'Select duration' dropdown set to 'One week', an 'End date' dropdown set to 'Yesterday', and a 'Period' field showing '29/10/2020 - 04/11/2020'. Below this is a 'Week Filter' section with radio buttons for 'Weekdays' (selected), 'Weekends', 'Mon', 'Tue', 'Wed', 'Thu', 'Fri', 'Sat', and 'Sun'. The 'Select reports to print' section features a list of reports with checkboxes and date range selectors:

Report	From	To
<input checked="" type="checkbox"/> All		
<input checked="" type="checkbox"/> Dashboard	29-10-2020	04-11-2020
<input checked="" type="checkbox"/> Daily Summary	29-10-2020	04-11-2020
<input checked="" type="checkbox"/> Day by Day Overview	29-10-2020	04-11-2020
<input checked="" type="checkbox"/> Sensor Overlay	29-10-2020	04-11-2020
<input checked="" type="checkbox"/> Trend Analysis	29-10-2020	04-11-2020
<input checked="" type="checkbox"/> Alert	29-10-2020	04-11-2020


A 'Create Report' button is located at the bottom right of the form.

Messages, profile and password

Click the down arrow at the top right corner, and a drop-down menu will appear.

If you see a red dot, it means you have received a message.

9.1 Messages

Click on **Messages** in the drop-down menu to go to the Messages screen.

If a patient accepted your connection request, denied your connection request, or deleted you from her/his connections, you will receive a message.


Messages, profile and password

9.2 Profile

Click on **Profile** in the drop-down menu to go to the Profile screen.

Here you can edit your personal profile.

9.3 Change Password

Click on **Change Password** in the drop-down menu to go to the Change Password screen.

Here you can set a new password for your account.

Messages, profile and password

Change Password [→ Home](#)

Old password

New password

Re-enter password

[Forgot your password?](#)

Messages, profile and password

Glossary

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- Specijalni karakteri: ~!@#% ^ & * () _ - + = ` { } | [] \ : ; ' & # € £ ¢ ¤ ¥ ¦ § ¨ © ª « ¬ ® ¯ ° ± ² ³ ´ µ ¶ · ¸ ¹ º » ¼ ½ ¾ ¿