

	Android	iPhone	iPad	Apple Watch	Windows	Mac
Pacijent	<p>TouchCare Sistem: EasyTouch App Samostalni CGM: EasySense App Kompatibilni uređaji: Samsung: Galaxy S7 Edge, Galaxy S8, Galaxy S9, Galaxy S10, Galaxy A9, Galaxy A40, Galaxy J3, Galaxy J5, Galaxy Note 10 LG: Google Nexus 5 SonyXperia: Z5 Xiaomi: MI 5X, MI 6, MI 8, MIX 2, MIX 2S, Note 3 Huawei: P9, P10, P20, P30, P30 Pro, Honor 7C, Honor 8, Honor V10, Honor 20, Mate 20, Mate 30, Nova 3e Oppo: R9, Reno Z Vivo: X9, X27, Y97 Podržani operativni sistemi: Android 4.3 i noviji</p> <p>EasyPatch APP Kompatibilni uređaji: Samsung: Galaxy S7 Edge, Galaxy S8, Galaxy S10e, Galaxy A40, Galaxy A9, Galaxy J3, Galaxy J5, Galaxy Note 10 LG: Google Nexus 5 Xiaomi: MI 5X Huawei: P9, P10, P20, P30, P30 Pro, Honor 7C, Honor 8, Honor V10, Honor 20, Mate 20, Mate 30, Nova 3e Oppo: Reno Z Vivo: X9 Podržani operativni sistemi: Android 8.0 i noviji</p>	<p>TouchCare System: EasyTouch App Stand-alone CGM: EasySense App Kompatibilni uređaji: iPhone: 5C, 5S, 6, 6 Plus, 6S, 6S Plus, SE, 7, 7 Plus, 8, X, XR, XS, 11 iPod touch: 5th Gen Podržani operativni sistemi: IOS 8.0 i noviji</p> <p>EasyPatch APP Kompatibilni uređaji: iPhone: 5C, 5S, 6, 6S, 7Plus, 8, X, XR, XS, 11 Podržani operativni sistemi: IOS 9.3.5 i noviji</p>	<p>TouchCare System: EasyTouch App Stand-alone CGM: EasySense App Kompatibilni uređaji: iPad: 3, 4, Air, Air 2, Mini, Mini 2, Mini 3, Mini 4, Pro</p>	<p>EasyTouch App EasySense App, EasyFollow App Kompatibilni uređaji: iPhone: 5S, 6, 6 Plus, 6S, 6S Plus, SE, 7, 7 plus, 8, XR Apple watch: series 1, 2, 3. Podržani operativni sistemi: IOS 11.0 i noviji, Watch OS 4.0 i noviji</p>	<p>EasyView personal website https://easyview.medtrum.eu OS: Windows 7, 8, 10 Pretraživač: Google Chrome 61.0 i noviji, Firefox 56.0.2 i noviji, IE 11 i noviji, Microsoft Edge 41 i noviji</p>	<p>EasyView personal website https://easyview.medtrum.eu OS: Mac OS 10.10.5 i noviji Pretraživač: Safari 8.0.8 i noviji</p>

Negovatelj	EasyFollow App Kompatibilni uređaji: Samsung: Galaxy S7 Edge, Galaxy S8, Galaxy S9, Galaxy S10, Galaxy A9, Galaxy A40, Galaxy J3, Galaxy J5, Galaxy Note10 LG: Google Nexus 5 SonyXperia: Z5 Xiaomi: MI 5X, MI 6, MI 8, MIX 2, MIX 2S, Note3 Huawei: P9, P10, P20, P30, P30 Pro, Honor 7C, Honor 8, Honor V10, Honor 20, Mate 20, Mate 30, Nova 3e Oppo: R9, Reno Z Vivo: X9, X27, Y97 Podržani operativni sistemi: Android 4.3 i noviji	EasyFollow App Kompatibilni uređaji: iPhone: 5C, 5S, 6, 6 Plus, 6S, 6S Plus, SE, 7, 7 plus, 8, X, XR, XS, 11 iPod touch: 5th Gen Podržani operativni sistemi: IOS 8.0 i noviji	EasyFollow App Kompatibilni uređaji: iPad: 3, 4, Air, Air 2, Mini, Mini 2, Mini 3, Mini 4, Pro		EasyView personal website https://easyview.medtrum.eu OS: Windows 7, 8, 10 Pretraživač: Google Chrome 61.0 i noviji, Firefox 56.0.2 i noviji, IE 11 i noviji, Microsoft Edge 41 i noviji	EasyView personal website https://easyview.medtrum.eu OS: Mac OS 10.10.5 i noviji Pretraživač: Safari 8.0.8 i noviji

Lekar			EasyView iPad App Kompatibilni uređaji: iPad: 3, 4, Air, Air 2, Mini, Mini 2, Mini 3, Mini 4, Pro		EasyView Pro website www.easyviewpro.com OS: Windows 7, 8, 10 Pretraživač: Google Chrome 61.0 i noviji, Firefox 56.0.2 i noviji, IE 11 i noviji, Microsoft Edge 41 i noviji	EasyView Pro website www.easyviewpro.com OS: Mac OS 10.10.5 i noviji Pretraživač: Safari 8.0.8 i noviji
--------------	--	--	---	--	--	--